

Automated External Defibrillators

ZOLL®



When they need CPR,
you need ZOLL.

Confidence that comes only from **Real CPR Help.**

AED Plus

The best support when rescuing a victim of cardiac arrest

The ZOLL AED Plus gives even infrequent rescuers the confidence they need to help save lives. Simple audio prompts and illustrations reinforce every step in the resuscitation process.



The one-piece adult electrode pad (CPR-D-padz®) with Real CPR Help technology makes placement of electrodes faster and more accurate.

“Push Harder.” “Good Compressions.” “Press Shock Button.”


Choosing ZOLL® for your next automated external defibrillator (AED) gives your patients yet another reason to trust you. As the brand hospitals choose most for professional defibrillators, ZOLL has earned the respect of users at every skill level. That’s because our AEDs do more than provide the most advanced defibrillation technology.

The AED Plus® and the AED Pro® are the only automated external defibrillators to provide Real CPR Help® with real-time feedback on the depth and rate of chest compressions during cardiopulmonary resuscitation (CPR). In addition to audio prompts, both the AED Plus and AED Pro feature on-screen instructions, and a bar gauge that shows how deep to push in order to optimize compression depth. This is an important advance, given the American Heart Association’s 2010 Guidelines and their recommendation for chest compressions at least two inches deep.



AED Pro


*A rugged, advanced AED
for professional users*



Chest compression
depth gauge



The ZOLL AED Pro, whether used in a hospital or in the field, is designed to perform in any environment, and is used by all types of rescuers. With Real CPR Help and See-Thru CPR®, it provides the advanced functionality that professional rescuers and services require.





AED Plus

- Real CPR Help provides real-time feedback for depth and rate of chest compressions during CPR
- Supports the entire Chain of Survival, with easy-to-understand audio prompts and illustrations
- Has lowest total cost of ownership, once installed, because CPR-D-padz and batteries last five years, if not used
- Rugged design that resists dust and direct water spray (IP55 rating)
- Powered by 10 lithium batteries, available from retail stores everywhere



AED Pro

- Real CPR Help provides real-time feedback for depth and rate of chest compressions during CPR
- See-Thru CPR filters CPR artifact to display patient's underlying ECG rhythm
- Advanced capabilities, including 3-lead patient monitoring and manual override
- The only AED that has an IP55 rating for dust and water resistance and passes the 1.5-meter drop test
- Battery and electrode compatibility with ZOLL's line of professional defibrillators

ZOLL Medical Corporation, an Asahi Kasei Group company, develops and markets medical devices and software solutions that help advance emergency care and save lives, while increasing clinical and operational efficiencies. With products for defibrillation and monitoring, circulation and CPR feedback, data management, fluid resuscitation, and therapeutic temperature management, ZOLL provides a comprehensive set of technologies that help clinicians, EMS and fire professionals, and lay rescuers treat victims needing resuscitation and critical care. For more information, visit www.zoll.com.

The Asahi Kasei Group is a diversified group of companies led by holding company Asahi Kasei Corp., with operations in the chemicals & fibers, homes & constructions materials, electronics and health care business sectors. Its health care operations include devices and systems for critical care, dialysis, therapeutic apheresis, transfusion, and manufacture of biotherapeutics, as well as pharmaceuticals, diagnostic reagents, and nutritional products. With more than 25,000 employees around the world, the Asahi Kasei Group serves customers in more than 100 countries. For more information, visit www.asahi-kasei.co.jp/asahi/en/.

**ADVANCING
RESUSCITATION.
TODAY.®**

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-348-9011
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.